

D: Körpermasstabelle
 E: Measurement size guide
 F: Tableau de mesures anatomiques
 I: Tabelle delle misure



VENOSAN® 3000 Cotton
VENOSAN® 4000
VENOSAN® 5000
VENOSAN® 7000

Umfangmasse/Circumference measurements/
 Mesures de circonférence/Misure circonferenza (cm)

	SMALL	MEDIUM	LARGE	X-LARGE
cB	20–23	23–26	26–29	29–31
cC	29–39	33–42	36–45	39–47
cD	28–38	32–41	35–44	38–46
cG	48–63	54–67	58–72	63–74

Längenmasse/Length/Longueur/Lunghezza (cm)

ID	AD/ ADH kurz/short	35–40
ID	AD/ ADH lang/long	37–42
IG	AGH kurz/short	63–72
IG	AGH lang/long	68–78
IK	AT/ ATU Materna kurz/short	67–74
IK	AGG/ AT/ ATU Mat. lang/long	74–85

VENOSAN® Perflor (nur/only AD)

Umfangmasse/Circumference measurements/
 Mesures de Circonférence/Misure circonferenza (cm)

	II	III	IV	V	VI
cB	20–22	22–24	24–26	26–28	28–30
cC	31–36	34–39	37–42	40–45	43–48
cD	29–34	32–37	35–40	38–43	41–46

Längenmasse/Length/Longueur/Lunghezza (cm)

ID	AD	39–44
----	----	-------

VENOSAN® 2000 (nur/only AD)

Umfangmasse/Circumference measurements/
 Mesures de Circonférence/Misure circonferenza (cm)

	SMALL		MEDIUM		LARGE	
	NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS
cB	19–23	19–23	23–26	23–26	26–29	26–29
cC	30–36	35–40	35–40	39–45	39–45	44–49
cD	28–34	33–38	33–38	37–43	37–43	42–47

Längenmasse/Length/Longueur/Lunghezza (cm)

ID	AD	40–44
ID	AD kurz/short	36–40

